
# ENSURE HIGH QUALITY OF PRODUCTS

## **CORE-Injector**



No.	Product Name	Working Length(mm)	Needle			
INO.			Gauze	O.D.(mm)	Thickness(mm)	Length(mm)
1	CORE Injector	2,300	23	0.6	0.1	4

- Smooth drug injection

**사**: 경상북도 경산시 와촌면 솔구불길 18,

천단메디컬융합섬유센터 501호 서울사무소 : 경기도 고양시 덕양구 통일로 140, 삼송테크노벨리 A동 222호

**Tel**. 02-866-3514 **Fax**. 02-6919-1346


INCORE Co., Ltd.

**Head office**: 501, 18, Solgubul-gil, Wachon-myeon, Gyeongsan-si, Gyeongsangbuk-do **Seoul office**: A-222, 140, Tongil-ro, Deogyang-gu, Goyang-si,

Gyeonggi-do

**Tel**. +82-2-866-3514 **Fax**. +82-2-6919-1346

Destributed By.





INCORE

INNOVATIVE CORE TECHNOLOGLES OF

## CORE-Forcep 위·대장 일회용 내시경 포셉

Features	Benefits		
Various Lineup	- Various cup sizes, open length 5.0mm / 6.0mm / 7.0mm		
Shape of Cup  Thinness of Cup →  Larger Specimen Volume  Double Cut Cup Edge	- Possible to withdraw more specimen compared to the same size - Perfect incision, Minimal bleeding		
<b>Double Wire</b> Stronger Force Delicate Operation	- Powerful force transmission to the cup while operating - Swing head, convenient operation		

No.	Product Name	Working	Сир			
INO.	Product Name	Length(mm)	Tooth Type	Open length		
1	CORE Forcep	1,600 (Gastro)	w/o Tooth	5.0 / 6.0 / 7.0		
2	CORE Forcep	2,300 (Colon)	w/o Tooth	5.0 / 6.0 / 7.0		

INNOVATIVE CORE TECHNOLOGLES OF MEDICAL DEVLCES



### **CORE-Snare**

위 · 대장 일회용 내시경 스네어



Features	Benefits		
Various Lineup	- Thickness of Loop / Softness of Loop / Multi Shape of Loop		
Advantage of Handle Ergonomic handle	- Prevention of slipperiness		
Advantage of Loop Delicate Operation	- Possible to predict the accurate open and close		

	Product Name	Loop			Sheath	
No.		Shape	Width(mm)	Strength	O.D.(mm)	Working Length(mm)
1	CORE Snare	Oval	10	Soft / Stiff	2.4	2,300
2	CORE Snare		15	Soft / Stiff	2.4	2,300
3	CORE Snare		24	Soft / Stiff	2.4	2,300
4	CORE Snare		30	Soft / Stiff	2.4	2,300
5	CORE Snare		10	Stiff	2.4	2,300
6	CORE Snare		15	Stiff	2.4	2,300
7	CORE Snare	Crescent	24	Stiff	2.4	2,300

