



# INCORE

Disposable Endoscopic Instruments

**INCORE Co., Ltd.**

Head Office 11, Hyeoksin-daero 78-gil, Dong-gu, Daegu, Republic of Korea

Branch Office (501), 18, Solgubul-gil, Wachon-myeon, Gyeongsan-si, Gyeongsangbuk-do, Republic of Korea

Tel. +82-2-866-3514 Fax. +82-2-6919-1346



# CORE-Forcep

01



Features	Benefits
Various Lineup	• Various cup size, Open length 5.0mm / 6.0mm / 7.0mm
Shape of Cup	• Perfect incision, Minimal bleeding
Double Wire	• Powerful force transmission to the cup while operating • Swing head, convenient operation

No.	Product Name	Working Length(mm)	Cup	
			Open Length(mm)	Size(mm)
1	CORE Forcep	1,600(Gastro)	5.0/6.0/7.0	1.7/2.0/2.3
2	CORE Forcep	2,300(Colon)	5.0/6.0/7.0	1.7/2.0/2.3

02

# CORE-Snare



Features	Benefits
Various Lineup	• Thickness of Loop/Softness of Loop/Multi Shape of Loop
Advantage of Handle	• Ergonomic handle • Prevention of slipperiness
Advantage of Loop	• Delicate Operation • Possible to predict the accurate open and close

No.	Type	Loop	Width(mm)	Sheath	
		Shape		O.D.(mm)	Working Length
1	Soft N/R	 Oval	10/15/24/30	2.4	2,300
2	Soft R			2.4	2,300
3	Stiff N/R			2.4	2,300
4	Stiff R			2.4	2,300
5	Hot/Cold N/R	 Diamond	10/15	2.4	2,300
6	Hot/Cold R			2.4	2,300

# CORE-Injector

03



Features	Benefits
One-Click Type	<ul style="list-style-type: none"> <li>• Possible to predict the accurate open and close</li> </ul>
Advantage of Outer Tube	<ul style="list-style-type: none"> <li>• No channel damage with smooth insertion</li> <li>• Prevention for projection of Inner Tube</li> <li>• Accurate reaching to the drug injection port</li> </ul>
Advantage of Needle	<ul style="list-style-type: none"> <li>• Delicate injection without slip</li> <li>• Smooth drug injection</li> </ul>

No.	Product Name	Working Length(mm)	Gauge	O.D.(mm)	Needle	
					Thickness(mm)	Length(mm)
1	CORE Injector	2,300	23	0.6	0.1	4
2	CORE Injector	2,300	25	0.6	0.1	4
3	CORE Injector	1,600	23	0.6	0.1	4
4	CORE Injector	1,600	25	0.6	0.1	4

04

# CORE-Clip




Features	Benefits
Rapid hemostatic effect by grasping opening site	<ul style="list-style-type: none"> <li>• It can be used immediately without connection with fixing device</li> </ul>
High clipping pressure	<ul style="list-style-type: none"> <li>• Clipping force is strong and enough to stop the bleeding</li> </ul>
Convenient function	<ul style="list-style-type: none"> <li>• Re-openable</li> <li>• Rotatable clipping head</li> </ul>

No.	Product Name	Sheath		Clip
		O.D.(mm)	W/Length	Open Length
1	Core-Clip	2.1	1,600	11
2	Core-Clip	2.1	1,600	14
3	Core-Clip	2.1	2,300	11
4	Core-Clip	2.1	2,300	14

# CORE-Knife

05



Features	Benefits
Various Types – D/I/T/L Types	Wide selection of users with various types
Controllable incision depth	Ensure safety during incision at constant drive intervals
Incision with less bleeding	Ensure safety against heat and bleeding by ceramic tip, the upper insulator

No.	Product Name	Sheath		Knife	
		O.D.	W/Length	Tip O.D.	Wire O.D.
1	CORE Knife D Type	2.3	2,300	2.00	0.43
2	CORE Knife D Type	2.3	1,800	2.00	0.43
3	CORE Knife I Type	2.3	2,300	1.00	0.43
4	CORE Knife I Type	2.3	1,800	1.00	0.43
5	CORE Knife T Type	2.3	2,300	0.60	0.43
6	CORE Knife T Type	2.3	1,800	0.60	0.43
7	CORE Knife L Type	2.3	2,300	1.65	0.43
8	CORE Knife L Type	2.3	1,800	1.65	0.43

06

# CORE-Hotbiopsy



Features	Benefits
Various Lineup	<ul style="list-style-type: none"> <li>• Various cup size, open length 5.0mm/7.0mm</li> <li>• Two types for Gastro and Colon</li> </ul>
Shape of Cup	Perfect incision, Minimal bleeding
Minimum electrical resistance	Excellent electrical performance makes excellent cutting performance.

No.	Product Name	Working Length(mm)	Cup	
			Open length(mm)	Size(mm)
1	CORE Hotbiopsy	1,600(Gastro)	5.0/7.0	1.7/2.3
2	CORE Hotbiopsy	2,300(Colon)	5.0/7.0	1.7/2.3